


MICROTEC LED FLUORESCENCE UNITS




The MICROTEC LED incident fluorescence unit is available in four forms and is designed to be used with the MICROTEC L3MED microscope. The following guide enables the appropriate unit to be selected.




This unit is recommended for use with the following fluorochromes:

- **DAPI**
- **HOECHST 33342/33258**
- **AMCA/AMCA-X**
- **Alexa Fluor 350**
- **Fluorescent Brightener Uvitex2B**
- **Calcofluor white**




This unit is recommended for use with the following fluorochromes:

- **Auramine O**
- **wtGFP**



This unit is recommended for use with the following fluorochromes:

- **EGFP**
- **FITC**
- **Alexa Fluor 488**
- **Cy 2**
- **DIO**
- **Fluo 4**
- **Acridine Orange**



This unit is recommended for use with the following fluorochromes:

- **Texas Red**
- **Texas Red-X**
- **Cy3.5**
- **Mito Tracker Red**
- **TRITC**